

Datasheet

Rabbit Anti-Chicken IgG F(ab')₂ secondary antibody

Catalog Number: PAB9178

Regulatory Status: For research use only (RUO)

Product Description: Rabbit anti-chicken IgG recognizes chicken IgG F(ab')₂ fragment. This secondary antibody was purified using delipidation and defibrination.

Immunogen: Chicken IgG F(ab')₂ fragment

Host: Rabbit

Reactivity: Chicken

Applications: ELISA, IHC, WB
(See our web site product page for detailed applications information)

Protocols: See our web site at <http://www.abnova.com/support/protocols.asp> or product page for detailed protocols

Specificity: No reaction was observed against chicken IgG F(c).

Form: Lyophilized

Conjugation: Unconjugated

Purification: This product was prepared from monospecific antiserum by a delipidation and defibrination. Assay by immunoelectrophoresis resulted in a single precipitin arc against chicken IgG F(ab')₂ and chicken serum.

Recommend Usage: The optimal working dilution should be determined by the end user.

Storage Buffer: Lyophilized from 20 mM KH₂PO₄, 150 mM NaCl, pH 7.2 (0.01% sodium azide)

Storage Instruction: Store at 4°C on dry atmosphere. After reconstitution with 2 mL of deionized water, store at -20°C or lower. Aliquot to avoid repeated freezing and thawing.