

Datasheet

Rabbit Anti-Dog IgG F(ab')₂ secondary antibody

Catalog Number: PAB9182

Regulatory Status: For research use only (RUO)

Product Description: Rabbit anti-dog IgG recognizes dog IgG F(ab')₂ fragment. This secondary antibody was purified using antigen affinity chromatography.

Immunogen: Dog IgG F(ab')₂ fragment

Host: Rabbit

Reactivity: Dog

Applications: ELISA, IHC, WB

(See our web site product page for detailed applications information)

Protocols: See our web site at

<http://www.abnova.com/support/protocols.asp> or product page for detailed protocols

Specificity: No reaction was observed against dog IgG F(ab')₂.

Form: Liquid

Conjugation: Unconjugated

Purification: This product was prepared from monospecific antiserum by immunoaffinity chromatography, followed by solid phase adsorption(s) to remove any unwanted reactivities.

Recommend Usage: The optimal working dilution should be determined by the end user.

Storage Buffer: In 20 mM KH₂PO₄, 150 mM NaCl, pH 7.2

Storage Instruction: Store at 4°C. For long term storage store at -20°C.
Aliquot to avoid repeated freezing and thawing.