

## Datasheet

### Mouse Anti-Goat IgG (H&L) secondary antibody

**Catalog Number:** PAB9210

**Regulatory Status:** For research use only (RUO)

**Product Description:** Mouse anti-Goat IgG recognizes Goat IgG whole molecule. This secondary antibody was purified using delipidation and defibrination.

**Immunogen:** Goat IgG whole molecule

**Host:** Mouse

**Reactivity:** Goat

**Applications:** ELISA, IHC, WB

(See our web site product page for detailed applications information)

**Protocols:** See our web site at

<http://www.abnova.com/support/protocols.asp> or product page for detailed protocols

**Form:** Lyophilized

**Conjugation:** Unconjugated

**Purification:** This product was prepared from monospecific antiserum by a delipidation and defibrination. Assay by immunoelectrophoresis resulted in a single precipitin arc against goat IgG and goat serum.

**Recommend Usage:** The optimal working dilution should be determined by the end user.

**Storage Buffer:** Lyophilized from 10 mM  $\text{KH}_2\text{PO}_4$ , 150 mM NaCl, pH 7.2

**Storage Instruction:** Store at 4°C on dry atmosphere. After reconstitution with 2 mL of deionized water, store at -20°C or lower. Aliquot to avoid repeated freezing and thawing.