

## Datasheet

### Rabbit Anti-Goat IgG F(c) secondary antibody

**Catalog Number:** PAB9227

**Regulatory Status:** For research use only (RUO)

**Product Description:** Rabbit anti-goat IgG recognizes goat IgG F(c) fragment. This secondary antibody was purified using delipidation and defibrination.

**Immunogen:** Goat IgG F(c) fragment

**Host:** Rabbit

**Reactivity:** Goat

**Applications:** ELISA, IHC, WB

(See our web site product page for detailed applications information)

**Protocols:** See our web site at

<http://www.abnova.com/support/protocols.asp> or product page for detailed protocols

**Specificity:** No reaction was observed against goat IgG F(ab')<sub>2</sub>.

**Form:** Lyophilized

**Conjugation:** Unconjugated

**Purification:** This product was prepared from monospecific antiserum by a delipidation and defibrination. Assay by immunoelectrophoresis resulted in a single precipitin arc against goat IgG F(c) and goat serum.

**Conjugation Note:**

**Recommend Usage:** The optimal working dilution should be determined by the end user.

**Storage Buffer:** Lyophilized from antiserum

**Storage Instruction:** Store at 4°C on dry atmosphere. After reconstitution with 2 mL of deionized water, store at -20°C or lower. Aliquot to avoid repeated freezing and thawing.