

Datasheet

Goat Anti-Golden Syrian hamster IgG (H&L) secondary antibody

Catalog Number: PAB9229

Regulatory Status: For research use only (RUO)

Product Description: Goat anti-Golden Syrian hamster IgG recognizes Golden Syrian hamster IgG whole molecule. This secondary antibody was purified using antigen affinity chromatography.

Immunogen: Golden Syrian hamster IgG whole molecule

Host: Goat

Reactivity: Hamster

Applications: ELISA, IHC, WB

(See our web site product page for detailed applications information)

Protocols: See our web site at <http://www.abnova.com/support/protocols.asp> or product page for detailed protocols

Specificity: No reaction was observed against bovine, chicken, goat, Guinea pig, horse, human, mouse, rabbit, rat and sheep serum proteins. Diminished reactivity will occur with Armenian hamster IgG.

Form: Liquid

Conjugation: Unconjugated

Purification: This product was prepared from monospecific antiserum by immunoaffinity chromatography, followed by solid phase adsorption(s) to remove any unwanted reactivities.

Recommend Usage: The optimal working dilution should be determined by the end user.

Storage Buffer: In 20 mM KH_2PO_4 , 150 mM NaCl, pH 7.2 (0.01% sodium azide)

Storage Instruction: Store at 4°C. For long term storage store at -20°C.
Aliquot to avoid repeated freezing and thawing.