

Datasheet

Goat Anti-Rabbit IgG (H&L) secondary antibody

Catalog Number: PAB9393

Regulatory Status: For research use only (RUO)

Product Description: Goat anti-rabbit IgG recognizes rabbit IgG whole molecule. This secondary antibody was purified using antigen affinity chromatography.

Immunogen: Rabbit IgG whole molecule

Host: Goat

Reactivity: Rabbit

Applications: ELISA, IHC, WB

(See our web site product page for detailed applications information)

Protocols: See our web site at

<http://www.abnova.com/support/protocols.asp> or product page for detailed protocols

Specificity: No reaction was observed against anti-Pepsin, anti-goat IgG F(c) or bovine, chicken, goat, Guinea pig, hamster, horse, human, mouse, rat and sheep serum proteins.

Form: Liquid

Conjugation: Unconjugated

Purification: This product was prepared from monospecific antiserum by immunoaffinity chromatography, followed by solid phase adsorption(s) to remove any unwanted reactivities, pepsin digestion and chromatographic separation.

Recommend Usage: The optimal working dilution should be determined by the end user.

Storage Buffer: In 20 mM KH₂PO₄, 150 mM NaCl, pH 7.2 (0.01% sodium azide)

Storage Instruction: Store at 4°C. For long term storage store at -20°C.

Aliquot to avoid repeated freezing and thawing.