

## Datasheet

### Goat Anti-Rabbit IgG (H&L) secondary antibody

**Catalog Number:** PAB9396

**Regulatory Status:** For research use only (RUO)

**Product Description:** Goat anti-rabbit IgG recognizes rabbit IgG whole molecule. This secondary antibody was purified using antigen affinity chromatography.

**Immunogen:** Rabbit IgG whole molecule

**Host:** Goat

**Reactivity:** Rabbit

**Applications:** ELISA, IHC, WB  
(See our web site product page for detailed applications information)

**Protocols:** See our web site at  
<http://www.abnova.com/support/protocols.asp> or product page for detailed protocols

**Specificity:** No reaction was observed against anti-Papain or anti-goat IgG F(c).

**Form:** Liquid

**Conjugation:** Unconjugated

**Purification:** This product was prepared from monospecific antiserum by immunoaffinity chromatography, followed by solid phase adsorption(s) to remove any unwanted reactivities, papain digestion and chromatographic separation.

**Recommend Usage:** The optimal working dilution should be determined by the end user.

**Storage Buffer:** In 20 mM  $\text{KH}_2\text{PO}_4$ , 150 mM NaCl, pH 7.2 (0.01% sodium azide)

**Storage Instruction:** Store at 4°C. Do not Freeze.