

Datasheet

Goat Anti-Rat IgG F(c) secondary antibody

Catalog Number: PAB9419

Regulatory Status: For research use only (RUO)

Product Description: Goat anti-rat IgG recognizes rat IgG F(c) fragment. This secondary antibody was purified using antigen affinity chromatography.

Immunogen: Rat IgG F(c) fragment

Host: Goat

Reactivity: Rat

Applications: ELISA, IHC, WB

(See our web site product page for detailed applications information)

Protocols: See our web site at

<http://www.abnova.com/support/protocols.asp> or product page for detailed protocols

Specificity: No reaction was observed against rat IgG F(ab')₂

Form: Liquid

Conjugation: Unconjugated

Purification: This product was prepared from monospecific antiserum by immunoaffinity chromatography, followed by solid phase adsorption(s) to remove any unwanted reactivities.

Recommend Usage: The optimal working dilution should be determined by the end user.

Storage Buffer: In 20 mM KH₂PO₄, 150 mM NaCl, pH 7.2 (0.01% sodium azide)

Storage Instruction: Store at 4°C. For long term storage store at -20°C.
Aliquot to avoid repeated freezing and thawing.