

## Datasheet

### Goat Anti-Rat IgG IgA IgM (H&L) secondary antibody

**Catalog Number:** PAB9420

**Regulatory Status:** For research use only (RUO)

**Product Description:** Goat anti-rat IgG recognizes rat IgG IgA IgM whole molecules. This secondary antibody was purified using antigen affinity chromatography.

**Immunogen:** Rat IgG, IgA and IgM whole molecules

**Host:** Goat

**Reactivity:** Rat

**Applications:** ELISA, IHC, WB  
(See our web site product page for detailed applications information)

**Protocols:** See our web site at <http://www.abnova.com/support/protocols.asp> or product page for detailed protocols

**Form:** Liquid

**Conjugation:** Unconjugated

**Purification:** This product was prepared from polyspecific antiserum by immunoaffinity chromatography, followed by solid phase adsorption(s) to remove any unwanted reactivities.

**Recommend Usage:** The optimal working dilution should be determined by the end user.

**Storage Buffer:** In 100 mM  $\text{Na}_2\text{B}_4\text{O}_7$ , 75 mM NaCl, 5 mM EDTA, pH 8.0 (0.01% sodium azide)

**Storage Instruction:** Store at 4°C. For long term storage store at -20°C.  
Aliquot to avoid repeated freezing and thawing.