

## Datasheet

### Goat Anti-Rat IgG (H&L) secondary antibody

**Catalog Number:** PAB9424

**Regulatory Status:** For research use only (RUO)

**Product Description:** Goat anti-rat IgG recognizes rat IgG whole molecule. This secondary antibody was purified using antigen affinity chromatography.

**Immunogen:** Rat IgG whole molecule

**Host:** Goat

**Reactivity:** Rat

**Applications:** ELISA, IHC, WB

(See our web site product page for detailed applications information)

**Protocols:** See our web site at

<http://www.abnova.com/support/protocols.asp> or product page for detailed protocols

**Specificity:** No reaction was observed against anti-Pepsin, anti-goat IgG F(c) or bovine, horse, human, mouse, rabbit and sheep serum proteins.

**Form:** Liquid

**Conjugation:** Unconjugated

**Purification:** This product was prepared from monospecific antiserum by immunoaffinity chromatography, followed by solid phase adsorption(s) to remove any unwanted reactivities, pepsin digestion and chromatographic separation.

**Recommend Usage:** The optimal working dilution should be determined by the end user.

**Storage Buffer:** In 20 mM  $\text{KH}_2\text{PO}_4$ , 150 mM NaCl, pH 7.2 (0.01% sodium azide)

**Storage Instruction:** Store at 4°C. For long term storage store at -20°C.  
Aliquot to avoid repeated freezing and thawing.