

## Datasheet

### Rabbit Anti-Rat IgG (H&L) secondary antibody

**Catalog Number:** PAB9435

**Regulatory Status:** For research use only (RUO)

**Product Description:** Rabbit anti-rat IgG recognizes rat IgG whole molecule. This secondary antibody was purified using antigen affinity chromatography.

**Immunogen:** Rat IgG whole molecule

**Host:** Rabbit

**Reactivity:** Rat

**Applications:** ELISA, IHC, WB  
(See our web site product page for detailed applications information)

**Protocols:** See our web site at <http://www.abnova.com/support/protocols.asp> or product page for detailed protocols

**Specificity:** No reaction was observed against anti-Pepsin and anti-rabbit IgG F(c).

**Form:** Liquid

**Conjugation:** Unconjugated

**Purification:** This product was prepared from monospecific antiserum by immunoaffinity chromatography, followed by solid phase adsorption(s) to remove any unwanted reactivities, pepsin digestion and chromatographic separation.

**Recommend Usage:** The optimal working dilution should be determined by the end user.

**Storage Buffer:** In 20 mM KH<sub>2</sub>PO<sub>4</sub>, 150 mM NaCl, pH 7.2 (0.01% sodium azide).

**Storage Instruction:** Store at 4°C. For long term storage store at -20°C.  
Aliquot to avoid repeated freezing and thawing.