

Datasheet

Rabbit Anti-Rat IgG (H&L) secondary antibody

Catalog Number: PAB9436

Regulatory Status: For research use only (RUO)

Product Description: Rabbit anti-rat IgG recognizes rat IgG whole molecule. This secondary antibody was purified using antigen affinity chromatography.

Immunogen: Rat IgG whole molecule

Host: Rabbit

Reactivity: Rat

Applications: ELISA, IHC, WB

(See our web site product page for detailed applications information)

Protocols: See our web site at

<http://www.abnova.com/support/protocols.asp> or product page for detailed protocols

Specificity: No reaction was observed against anti-Pepsin, anti-rabbit IgG F(c) or bovine, horse, human, mouse, rabbit & sheep serum proteins.

Form: Liquid

Conjugation: Unconjugated

Purification: This product was prepared from monospecific antiserum by immunoaffinity chromatography, followed by solid phase adsorption(s) to remove any unwanted reactivities, pepsin digestion and chromatographic separation.

Recommend Usage: The optimal working dilution should be determined by the end user.

Storage Buffer: In 20 mM KH_2PO_4 , 150 mM NaCl, pH 7.2 (0.01% sodium azide)

Storage Instruction: Store at 4°C. For long term storage store at -20°C.
Aliquot to avoid repeated freezing and thawing.