

Datasheet

Rabbit Anti-Rat IgG (H&L) secondary antibody

Catalog Number: PAB9441

Regulatory Status: For research use only (RUO)

Product Description: Rabbit anti-rat IgG recognizes rat IgG whole molecule. This secondary antibody was purified using delipidation and defibrination.

Immunogen: Rat IgG whole molecule

Host: Rabbit

Reactivity: Rat

Applications: ELISA, IHC, WB
(See our web site product page for detailed applications information)

Protocols: See our web site at <http://www.abnova.com/support/protocols.asp> or product page for detailed protocols

Form: Lyophilized

Conjugation: Unconjugated

Purification: This product was prepared from monospecific antiserum by a delipidation and defibrination. Assay by immunoelectrophoresis resulted in a single precipitin arc against rat IgG and rat serum.

Conjugation Note:

Recommend Usage: The optimal working dilution should be determined by the end user.

Storage Buffer: Lyophilized from 10 mM sodium phosphate, 140 mM NaCl, pH 7.4

Storage Instruction: Store at 4°C on dry atmosphere. After reconstitution with 2 mL of deionized water, store at -20°C or lower. Aliquot to avoid repeated freezing and thawing.