

Datasheet

Rabbit Anti-Pig IgG F(ab')₂ secondary antibody

Catalog Number: PAB9467

Regulatory Status: For research use only (RUO)

Product Description: Rabbit anti-pig IgG recognizes pig IgG F(ab')₂ fragment. This secondary antibody was purified using antigen affinity chromatography.

Immunogen: Pig IgG F(ab')₂ fragment

Host: Rabbit

Reactivity: Pig

Applications: ELISA, IHC, IP, WB
(See our web site product page for detailed applications information)

Protocols: See our web site at <http://www.abnova.com/support/protocols.asp> or product page for detailed protocols

Form: Liquid

Conjugation: Unconjugated

Purification: This product was prepared from monospecific antiserum by immunoaffinity chromatography, followed by solid phase adsorption(s) to remove any unwanted reactivities.

Recommend Usage: The optimal working dilution should be determined by the end user.

Storage Buffer: In 20 mM KH₂PO₄, 150 mM NaCl, pH 7.2 (0.01% sodium azide)

Storage Instruction: Store at 4°C. For long term storage store at -20°C.
Aliquot to avoid repeated freezing and thawing.