

Datasheet

Rabbit Anti-Pig IgG F(ab')₂ secondary antibody

Catalog Number: PAB9474

Regulatory Status: For research use only (RUO)

Product Description: Rabbit anti-pig IgG recognizes pig IgG F(ab')₂ fragment. This secondary antibody was purified using delipidation and defibrination.

Immunogen: Pig IgG F(ab')₂ fragment

Host: Rabbit

Reactivity: Pig

Applications: ELISA, IHC, WB

(See our web site product page for detailed applications information)

Protocols: See our web site at

<http://www.abnova.com/support/protocols.asp> or product page for detailed protocols

Form: Lyophilized

Conjugation: Unconjugated

Purification: This product was prepared from monospecific antiserum by a delipidation and defibrination. Assay by immunoelectrophoresis resulted in a single precipitin arc against pig IgG F(ab')₂ and pig serum.

Conjugation Note:

Recommend Usage: The optimal working dilution should be determined by the end user.

Storage Buffer: Lyophilized from 10 mM sodium phosphate, 250 mM NaCl, pH 7.4 (0.01% (w/v) sodium azide)

Storage Instruction: Store at 4°C on dry atmosphere. After reconstitution with 2 mL of deionized water, store at -20°C or lower. Aliquot to avoid repeated freezing and thawing.