

## Datasheet

### BrdU polyclonal antibody

**Catalog Number:** PAB9791

**Regulatory Status:** For research use only (RUO)

**Product Description:** Sheep polyclonal antibody raised against BrdU.

**Immunogen:** BrdU conjugated with KLH.

**Host:** Sheep

**Applications:** IHC, IP, WB

(See our web site product page for detailed applications information)

**Protocols:** See our web site at

<http://www.abnova.com/support/protocols.asp> or product page for detailed protocols

**Form:** Liquid

**Recommend Usage:** Western Blot (25-100 ug/mL)

Immunoprecipitation (25-100 ug/mL)

Immunohistochemistry (25-100 ug/mL)

The optimal working dilution should be determined by the end user.

**Storage Buffer:** In PBS (0.08% sodium azide)

**Storage Instruction:** Store at -20°C.

Aliquot to avoid repeated freezing and thawing.