

## Datasheet

### Propidium Iodide

**Catalog Number:** U0332

**Regulation Status:** For research use only (RUO)

**Product Description:** Propidium iodide is a red-fluorescent nuclear and chromosome counterstain. Since propidium iodide is not permeant to live cells, it is also commonly used to detect dead cells in a population.

**Applications:** Cell Staining  
(See our web site product page for detailed applications information)

**Protocols:** See our web site at <http://www.abnova.com/support/protocols.asp> or product page for detailed protocols

**Storage Instruction:** Store at room temperature and protect from light.