

Monoclonal Antibody to Phosphotyrosine - PE -

Catalog No.:	AM00199RP-N
Quantity:	0.1 mg
Concentration:	0.1 mg/ml
Host / Isotype:	Mouse / IgG2b
Clone:	PY20
Format:	State: Liquid purified Ig Buffer System: PBS, containing 0.09% Sodium azide Label: PE – Conjugated to R. Phycoerythrin (RPE)
Applications:	Flow Cytometry. Other applications not tested. Optimal dilutions are dependent on conditions and should be determined by the user.
Specificity:	Reacts to Phosphotyrosine.
Storage:	Store the antibody undiluted at 2-8°C. Do not freeze! This product is photosensitive and should be protected from light. Shelf life: one year from despatch.