

Monoclonal Antibody to Mouse IgE (Epsilon heavy chain) - PE

Alternate names:	Mouse Immunoglobulin E
Catalog No.:	AM08013RP-S
Quantity:	0.1 mg
Concentration:	0.1 mg/ml
Host / Isotype:	Rat / IgG1
Clone:	23G3
Format:	State: Liquid purified Ig fraction. Purification: PBS containing 0.09% Sodium Azide as preservative and a stabilizing agent. Label: PE – R-phycoerythrin
Applications:	Flow Cytometry: < / = 0.2 µg/10e6 cells. Identification and enumeration of IgE+ cells. Immunofluorescence Microscopy: Identification and enumeration of IgE+ cells. Other applications not tested. Optimal dilutions are dependent on conditions and should be determined by the user.
Specificity:	This antibody reacts with the epsilon heavy chain of BALB/c mouse IgE as demonstrated by ELISA. Species: Mouse. Other species not tested.
Storage:	Store the antibody undiluted at 2-8°C. DO NOT FREEZE! This product is photosensitive and should be protected from light. Avoid repeated freezing and thawing. Shelf life: one year from despatch.