

Monoclonal Antibody to Mouse Kappa Light Chain - FITC

Catalog No.: AM08020FC-N

Quantity: 0.5 mg

Concentration: 0.5 mg/ml

Host / Isotype: Rat / IgG1

Clone: H139-52.1

Format: **State:** Liquid purified Ig fraction.
Buffer System: PBS containing 0.09% Sodium Azide as preservative.
Label: FITC – Fluorescein Isothiocyanate Isomer 1

Applications: **Flow Cytometry:** < / = 1 µg/10e6 cells.
Identification and enumeration of kappa+ cells.
Immunofluorescence Microscopy: Identification and enumeration of kappa+ cells.
Other applications not tested. Optimal dilutions are dependent on conditions and should be determined by the user.

Specificity: This antibody reacts with Mouse kappa light chains.
Species: Mouse.
Other species not tested.

Storage: Store the antibody undiluted at 2-8°C for one month or (in aliquots) at -20°C for longer.
This product is photosensitive and should be protected from light.
Avoid repeated freezing and thawing.
Shelf life: one year from despatch.