

Monoclonal Antibody to Rabbit IgG - FITC

Alternate names:	Rabbit Immunoglobulin G
Catalog No.:	AM08102FC-N
Quantity:	0.5 mg
Concentration:	0.5 mg/ml
Host / Isotype:	Mouse / IgG1
Clone:	2A9
Format:	State: Liquid purified Ig fraction. Buffer System: PBS containing 0.09% Sodium Azide as preservative. Label: FITC – Fluorescein Isothiocyanate Isomer 1
Applications:	Flow Cytometry: < / = 3 µg/10e6 cells. Other applications not tested. Optimal dilutions are dependent on conditions and should be determined by the user.
Specificity:	This antibody reacts with the heavy chain of all subclasses of rabbit IgG as demonstrated by ELISA. It may also react with IgG from other species. Species: Rabbit. Other species not tested.
Storage:	Store the antibody undiluted at 2-8°C for one month or in (aliquots) at -20°C for longer. This product is photosensitive and should be protected from light. Avoid repeated freezing and thawing. Shelf life: one year from despatch.

For research and in vitro use only. Not for diagnostic or therapeutic work.

Material Safety Datasheets are available at www.acris-antibodies.com or on request.

Antibody Hotline - Technical Questions - Antibody Location Service
Free Call: 0800-2274746 (Germany only) - www.acris-antibodies.com