

Monoclonal Antibody to Rabbit IgG - PE

Alternate names:	Rabbit Immunoglobulin G
Catalog No.:	AM08102RP-N
Quantity:	0.1 mg
Concentration:	0.1 mg/ml
Host / Isotype:	Mouse / IgG1
Clone:	2A9
Format:	State: Liquid purified Ig fraction. Buffer System: PBS containing 0.09% Sodium Azide as preservative and a stabilizing agent. Label: PE – R-Phycoerythrin
Applications:	Flow Cytometry: < / = 0.3 µg/10e6 cells. Other applications not tested. Optimal dilutions are dependent on conditions and should be determined by the user.
Specificity:	This antibody reacts with the heavy chain of all subclasses of rabbit IgG as demonstrated by ELISA. It may also react with IgG from other species. Species: Rabbit. Other species not tested.
Storage:	Store the antibody undiluted at 2-8°C. DO NOT FREEZE! This product is photosensitive and should be protected from light. Avoid repeated freezing and thawing. Shelf life: one year from despatch.