

Monoclonal Antibody to CD3 - PE

Alternate names:	T-cell surface antigen T3/Leu-4, T-cell surface glycoprotein CD3, T3/Leu-4
Catalog No.:	AM08105RP-N
Quantity:	0.1 mg
Concentration:	0.1 mg/ml
Background:	CD3 epsilon, a member of the immunoglobulin superfamily of cell surface receptors, is comprised of five invariable chains (designated as gamma, delta, epsilon, zeta and eta) ranging in size from 16-28 kDa. CD3 epsilon is closely associated with the T cell antigen receptor (TCR), and is present on thymocytes, CD8+ and CD4+ positive cells. The CD3 epsilon complex mediates signal transduction. It plays a major role in signaling during antigen recognition, leading to T cell activation. The TCR/CD3E complex of T lymphocytes consists of either a TCR alpha/beta or TCR gamma/delta heterodimer coexpressed at the cell surface with the invariant subunits of CD3 labeled gamma, delta, epsilon, zeta, and eta. CD3 epsilon is expressed on all T cells of all mouse strains.
Uniprot ID:	Q7YRN2
NCBI:	NP_999392.1
GeneID:	397455
Host / Isotype:	Mouse / IgG2b
Clone:	BB23-8E6
Format:	State: Liquid purified Ig fraction. Buffer System: PBS containing 0.09% Sodium Azide as preservative and a stabilizing agent. Label: PE – R-Phycoerythrin
Applications:	Flow Cytometry: < / = 0.2 µg/10e6 cells. (Ref.1) Other applications not tested. Optimal dilutions are dependent on conditions and should be determined by the user.
Specificity:	This antibody recognizes Porcine CD3e (Mr. 20 kDa). The Monoclonal antibody BB23-8E6 is Mitogenic when presented to peripheral blood mononuclear cells in immobilized form. (Ref.1) Species: Pig. Other species not tested.
Storage:	Store the antibody undiluted at 2-8°C. DO NOT FREEZE! This product is photosensitive and should be protected from light. Avoid repeated freezing and thawing. Shelf life: one year from despatch.

General Readings: 1. Pescovitz, M.D. Personal communication.

For research and in vitro use only. Not for diagnostic or therapeutic work.

Material Safety Datasheets are available at www.acris-antibodies.com or on request.

Antibody Hotline - Technical Questions - Antibody Location Service
Free Call: 0800-2274746 (Germany only) - www.acris-antibodies.com