

Polyclonal Antibody to Bromodeoxyuridine (BrdU) - Purified

Catalog No.: AP05008PU-N

Quantity: 0.25 mg

Concentration: Lot specific

Background: Bromodeoxyuridine was coupled to keyhole limpet hemocyanin (KLH) and was used to immunize sheep.

Host / Isotype: Sheep / IgG

Format: **State:** Liquid purified Ig (0.2µm sterile filtered)

Purification: Protein G purified

Buffer System: Phosphate buffered saline with 0.08% sodium azide

Applications: Western blotting: 25 - 100 µg/ml.
Immunoprecipitation: 25 - 100 µg/ml.
Immunohistochemistry: 25 - 100 µg/ml.
Other applications not tested. Optimal dilutions are dependent on conditions and should be determined by the user.

Specificity: This antibody reacts to BrdU.

Add. Information: Denature DNA sample first so that bases are accessible to anti BrdU antibody.

Storage: Store the antibody at -20°C.
Avoid repeated freezing and thawing.
Shelf life: one year from despatch.

Pictures:

