

Polyclonal Antibody to Poliovirus (Type 1, 2, 3) - FITC

Catalog No.:	AP05342FC-N
Quantity:	1 ml
Concentration:	4 - 5 mg/ml (P)
Host / Isotype:	Goat / IgG
Format:	State: Liquid Ig fraction Buffer System: Phosphate buffered saline containing 0.09% Sodium Azide (NaN ₃) and 1%w/v Bovine Serum Albumin Label: FITC – Fluorescein Isothiocyanate Isomer 1
Applications:	Immunofluorescence. Other applications not tested. Optimal dilutions are dependent on conditions and should be determined by the user.
Specificity:	This antibody reacts to Polivirus 1, 2 and 3.
Storage:	Store the antibody undiluted at 2-8°C for one month or (in aliquots) at -20°C for longer. Avoid repeated freezing and thawing. This product is photosensitive and should be protected from light. Shelf life: one year from despatch.
Caution:	(A full Health and Safety assessment is available upon request) This product contains Sodium Azide: a POISONOUS AND HAZARDOUS SUBSTANCE, which should be handled by trained staff only.