

Polyclonal Antibody to Hamster IgG (F(ab)2 specific) - FITC

Alternate names:	Hamster Immunoglobulin G
Catalog No.:	AP05374FC-N
Quantity:	0.4 mg
Concentration:	0.8 mg/ml
Host / Isotype:	Goat / IgG
Immunogen:	Hamster IgG
Format:	State: Liquid F(ab') ₂ fragmend Purification: Affinity chromatography Buffer System: Phosphate buffered saline containing 0.09% Sodium Azide, 0.5% Bovine Serum Albumin Label: FITC – Fluorescein Isothiocyanate Isomer 1
Applications:	Flow Cytometry: 1/100 - 1/200. Other applications not tested. Optimal dilutions are dependent on conditions and should be determined by the user.
Specificity:	This antibody recognises hamster IgG. The reagent has been adsorbed against both mouse and rat immunoglobulins to minimise cross-reactivity. Species: Hamster Other species not tested.
Storage:	Store the antibody undiluted at 2-8°C for one month or (in aliquots) at -20°C for longer. Avoid repeated freezing and thawing. This product is photosensitive and should be protected from light. Shelf life: one year from despatch.
Caution:	(A full Health and Safety assessment is available upon request) This product contains sodium azide: a POISONOUS AND HAZARDOUS SUBSTANCE which should be handled by trained staff only.