

Polyclonal Antibody to Rat IgG1 - FITC

Alternate names:	Rat Immunoglobulin G1
Catalog No.:	AP05378FC-N
Quantity:	0.5 mg
Concentration:	1.0 mg/ml
Host / Isotype:	Goat / IgG
Format:	State: Liquid purified IgG Purification: Affinity chromatography Buffer System: Phosphate buffered saline pH7.2 containing 0.09% Sodium Azide and 0.2% Bovine Serum Albumin Label: FITC – Fluorescein Isothiocyanate Isomer 1
Applications:	Flow Cytometry: 1/400. Other applications not tested. Optimal dilutions are dependent on conditions and should be determined by the user.
Specificity:	This antibody recognises IgG1. The antibody has been shown to react specifically with rat IgG1 by immunoelectrophoresis and ELISA. The antibody may cross react with IgG1 from other species. Species: Rat. Other species not tested.
Storage:	Store the antibody undiluted at 2-8°C for one month or (in aliquots) at -20°C for longer. Avoid repeated freezing and thawing. This product is photosensitive and should be protected from light. Shelf life: one year from despatch.

For research and in vitro use only. Not for diagnostic or therapeutic work.

Material Safety Datasheets are available at www.acris-antibodies.com or on request.

Antibody Hotline - Technical Questions - Antibody Location Service
Free Call: 0800-2274746 (Germany only) - www.acris-antibodies.com