

Polyclonal Antibody to Rat IgG1 - HRP

Alternate names:	Rat Immunoglobulin G1
Catalog No.:	AP05378HR-N
Quantity:	0.5 mg
Concentration:	1.0 mg/ml
Host / Isotype:	Goat / IgG
Format:	State: Liquid purified IgG Purification: Affinity chromatography Buffer System: Phosphate buffered saline pH7.2 containing 0.1% Proclin 300 and 0.2% Bovine Serum Albumin Label: HRP – Horseradish Peroxidase
Applications:	ELISA: 1:10,000. Western Blot: 1:1,000 - 10,000. Other applications not tested. Optimal dilutions are dependent on conditions and should be determined by the user.
Specificity:	This antibody recognises IgG1. The antibody has been shown to react specifically with rat IgG1 by immunoelectrophoresis and ELISA. The antibody may cross react with IgG1 from other species. Species: Rat. Other species not tested.
Storage:	Store the antibody undiluted at 2-8°C for one month or (in aliquots) at -20°C for longer. Avoid repeated freezing and thawing. Shelf life: one year from despatch.
Caution:	(A full Health and Safety assessment is available upon request) This product contains Sodium Azide: a POISONOUS AND HAZARDOUS SUBSTANCE, which should be handled by trained staff only.

For research and in vitro use only. Not for diagnostic or therapeutic work.

Material Safety Datasheets are available at www.acris-antibodies.com or on request.

Antibody Hotline - Technical Questions - Antibody Location Service
Free Call: 0800-2274746 (Germany only) - www.acris-antibodies.com