

Polyclonal Antibody to Rat IgG1 - Purified

Alternate names:	Rat Immunoglobulin G1
Catalog No.:	AP05378PU-N
Quantity:	0.5 mg
Concentration:	1.0 mg/ml
Host / Isotype:	Goat / IgG
Format:	State: Liquid purified IgG Purification: Affinity chromatography Buffer System: Phosphate buffered saline pH7.2 containing 0.09% Sodium Azide
Applications:	ELISA: 1µg/ml - 10µg/ml. Western Blot. Flow Cytometry. Other applications not tested. Optimal dilutions are dependent on conditions and should be determined by the user.
Specificity:	This antibody recognises IgG1. The antibody has been shown to react specifically with rat IgG1 by immunoelectrophoresis and ELISA. The antibody may cross react with IgG1 from other species. Species: Rat. Other species not tested.
Storage:	Store the antibody undiluted at 2-8°C for one month or (in aliquots) at -20°C for longer. Avoid repeated freezing and thawing. Shelf life: one year from despatch.
Caution:	(A full Health and Safety assessment is available upon request) This product contains Sodium Azide: a POISONOUS AND HAZARDOUS SUBSTANCE, which should be handled by trained staff only.