

Polyclonal Antibody to Rat IgG2a - HRP

Alternate names:	Rat Immunoglobulin G2a
Catalog No.:	AP05379HR-N
Quantity:	0.5 mg
Concentration:	1.0 mg/ml
Host / Isotype:	Goat / IgG
Format:	State: Liquid purified IgG Purification: Affinity chromatography Buffer System: Phosphate buffered saline pH7.2 containing 0.1% ProClin300 and 0.2% Bovine Serum Albumin Label: HRP – Horseradish Peroxidase
Applications:	ELISA: 1:10,000. Western Blot: 1:1000 - 10000. Other applications not tested. Optimal dilutions are dependent on conditions and should be determined by the user.
Specificity:	This antibody recognises rat IgG2a. The antibody has been shown to react specifically with rat IgG2a by immunoelectrophoresis and ELISA. The antibody may cross react with IgG2a from other species. Species: Rat. Other species not tested.
Storage:	Store the antibody undiluted at 2-8°C for one month or (in aliquots) at -20°C for longer. Avoid repeated freezing and thawing. Shelf life: one year from despatch.
Caution:	(A full Health and Safety assessment is available upon request) This product contains Sodium Azide: a POISONOUS AND HAZARDOUS SUBSTANCE, which should be handled by trained staff only.