

Polyclonal Antibody to Mouse IgG1 - PE

Alternate names:	Mouse Immunoglobulin G1
Catalog No.:	AP05384RP-N
Quantity:	0.1 mg
Concentration:	0,5 mg/ml
Host / Isotype:	Goat / IgG
Immunogen:	Mouse IgG1 paraproteins
Format:	State: Liquid purified IgG Buffer System: Phosphate buffered saline pH7.4 containing 0,09% sodium azide and sucrose Label: PE – R.Phycoerythrin (RPE)
Applications:	Flow Cytometry: 1/20 - 1/30. Other applications not tested. Optimal dilutions are dependent on conditions and should be determined by the user.
Specificity:	This antibody recognises IgG1. This antibody has been cross absorbed against mouse IgM, IgG2a, IgG2b, IgG3 and IgA, pooled human sera and purified human paraproteins. The antibody shows minimal cross-reactivity with human immunoglobulins. Species: Mouse. Other species not tested.
Storage:	Store the antibody undiluted at 2-8°C. DO NOT FREEZE! This product is photosensitive and should be protected from light. Shelf life: 6 month from despatch.
Caution:	(A full Health and Safety assessment is available upon request) This product contains Thiomersal: a POISONOUS AND HAZARDOUS SUBSTANCE, which should be handled by trained staff only.