

Polyclonal Antibody to Mouse IgG2a - FITC

Alternate names:	Mouse Immunoglobulin 2a
Catalog No.:	AP05385FC-N
Quantity:	0.5 mg
Concentration:	1.0mg/ml
Host / Isotype:	Goat / IgG
Immunogen:	Mouse IgG2a paraproteins
Format:	State: Liquid purified IgG Buffer System: PBS pH7.4 containing 0,09% sodium azide Label: FITC – Fluorescein Isothiocyanate Isomer 1
Applications:	Flow Cytometry: 1/10. Immunofluorescence: 1/10. Other applications not tested. Optimal dilutions are dependent on conditions and should be determined by the user.
Specificity:	This antibody recognises Mouse IgG2a. This antibody has been cross absorbed against mouse IgM, IgG1, IgG2b, IgG3 and IgA, pooled human sera and purified human paraproteins. The antibody shows minimal cross-reactivity with human immunoglobulins. Species: Mouse. Other species not tested.
Storage:	Store the antibody undiluted at 2-8°C for one month or (in aliquots) at -20°C for longer. Avoid repeated freezing and thawing. This product is photosensitive and should be protected from light. Shelf life: one year from despatch.
Caution:	(A full Health and Safety assessment is available upon request) This product contains sodium azide: a POISONOUS AND HAZARDOUS SUBSTANCE which should be handled by trained staff only.