

Polyclonal Antibody to Mouse IgG2a - PE

Alternate names:	Mouse Immunoglobulin 2a
Catalog No.:	AP05385RP-N
Quantity:	0.1 mg
Concentration:	0,5 mg/ml
Host / Isotype:	Goat / IgG
Immunogen:	Mouse IgG2a paraproteins
Format:	State: Liquid purified IgG Buffer System: PBS, pH 7,4 containing 0,09% sodium azide and Sucrose Label: PE – R.Phycoerythrin (RPE)
Applications:	Flow Cytometry: 1/20 - 1/30. Other applications not tested. Optimal dilutions are dependent on conditions and should be determined by the user.
Specificity:	This antibody recognises Mouse IgG2a. This antibody has been cross absorbed against mouse IgM, IgG1, IgG2b, IgG3 and IgA, pooled human sera and purified human paraproteins. The antibody shows minimal cross-reactivity with human immunoglobulins. Species: Mouse. Other species not tested.
Storage:	Store the antibody undiluted at 2-8°C. DO NOT FREEZE! This product is photosensitive and should be protected from light. Shelf life: 6 month from despatch.
Caution:	(A full Health and Safety assessment is available upon request) This product contains sodium azide: a POISONOUS AND HAZARDOUS SUBSTANCE which should be handled by trained staff only.

For research and in vitro use only. Not for diagnostic or therapeutic work.

Material Safety Datasheets are available at www.acris-antibodies.com or on request.

Antibody Hotline - Technical Questions - Antibody Location Service
Free Call: 0800-2274746 (Germany only) - www.acris-antibodies.com