

Polyclonal Antibody to Human IgE - FITC

Alternate names:	Immunoglobulin E human
Catalog No.:	AP05393FC-N
Quantity:	1 mg
Concentration:	1.0 mg/ml
Host / Isotype:	Goat / IgG
Format:	State: Liquid purified IgG fraction. Purification: Affinity Chromatography. Buffer System: Phosphate buffered saline containing 0.09% Sodium Azide as preservative and 0.2% BSA as stabilizer. Label: FITC – Fluorescein Isothiocyanate Isomer 1
Applications:	Flow Cytometry. Other applications not tested. Optimal dilutions are dependent on conditions and should be determined by the user.
Specificity:	This antibody specifically recognises Human IgE. Species: Human and Rhesus Monkey. Other species not tested.
Storage:	Store the antibody undiluted at 2-8°C for one month or (in aliquots) at -20°C for longer. This product is photosensitive and should be protected from light. Avoid repeated freezing and thawing. Shelf life: one year from despatch.