

Polyclonal Antibody to Rabbit IgA - FITC

Alternate names:	Rabbit Immunoglobulin A
Catalog No.:	AP05550FC-N
Quantity:	1 mg
Concentration:	1.0 mg/ml
Host / Isotype:	Goat / IgG
Immunogen:	Purified rabbit IgA
Format:	State: Liquid purified IgG Buffer System: Phosphate buffered saline containing 0.09% Sodium Azide, 0.2% Bovine Serum Albumin Label: FITC – Fluorescein Isothiocyanate Isomer 1
Applications:	Immunohistochemistry on frozen sections: 1/200 - 1/2000. Flow Cytometry. Other applications not tested. Optimal dilutions are dependent on conditions and should be determined by the user.
Specificity:	This antibody recognises IgA and shows no cross - reactivity with other rabbit immunoglobulin classes in immunoelectrophoresis. Species: Rabbit. Other species not tested.
Storage:	Store the antibody undiluted at 2-8°C for one month or (in aliquots) at -20°C for longer. Avoid repeated freezing and thawing. This product is photosensitive and should be protected from light. Shelf life: one year from despatch.
Caution:	(A full Health and Safety assessment is available upon request) This product contains sodium azide: a POISONOUS AND HAZARDOUS SUBSTANCE which should be handled by trained staff only.

For research and in vitro use only. Not for diagnostic or therapeutic work.

Material Safety Datasheets are available at www.acris-antibodies.com or on request.

Antibody Hotline - Technical Questions - Antibody Location Service
Free Call: 0800-2274746 (Germany only) - www.acris-antibodies.com