

## Polyclonal Antibody to Turkey IgG (H+L chain) - FITC

<b>Alternate names:</b>	Turkey Immunoglobulin G
<b>Catalog No.:</b>	AP08138FC-N
<b>Quantity:</b>	1 mg
<b>Concentration:</b>	1.0 mg/ml
<b>Host / Isotype:</b>	Goat / IgG
<b>Format:</b>	<b>State:</b> Liquid purified IgG fraction. <b>Purification:</b> Affinity Chromatography on Chicken IgG covalently linked to agarose. <b>Buffer System:</b> PBS containing 0.09% Sodium Azide as preservative. <b>Label:</b> FITC – Fluorescein Isothiocyanate Isomer 1
<b>Applications:</b>	<b>Immunofluorescence:</b> < / = 1 µg/10e6 cells. <b>ELISA.</b> <b>Western blotting.</b> <b>Dot- and Slot-Immunoblotting.</b> <b>Immunohistochemistry.</b> Other applications not tested. Optimal dilutions are dependent on conditions and should be determined by the user.
<b>Specificity:</b>	Reacts with the heavy and light chains of turkey IgG as demonstrated by ELISA. May also react with the light chains of other turkey immunoglobulins. <b>Species:</b> Turkey. Other species not tested.
<b>Storage:</b>	Store the antibody undiluted at 2-8°C for one month or in (aliquots) at -20°C for longer. This product is photosensitive and should be protected from light. Avoid repeated freezing and thawing. Shelf life: one year from despatch.

**For research and in vitro use only. Not for diagnostic or therapeutic work.**

Material Safety Datasheets are available at [www.acris-antibodies.com](http://www.acris-antibodies.com) or on request.

Antibody Hotline - Technical Questions - Antibody Location Service  
Free Call: 0800-2274746 (Germany only) - [www.acris-antibodies.com](http://www.acris-antibodies.com)