

Polyclonal Antibody to Turkey IgG (H+L chain) - PE

Alternate names:	Turkey Immunoglobulin G
Catalog No.:	AP08138RP-N
Quantity:	0.5 mg
Concentration:	0.5 mg/ml
Host / Isotype:	Goat / IgG
Format:	State: Liquid purified IgG fraction. Purification: Affinity Chromatography on Chicken IgG covalently linked to agarose. Buffer System: PBS containing 0.09% Sodium Azide as preservative and a stabilizing agent. Label: PE – R-Phycoerythrin
Applications:	Immunofluorescence. ELISA. Western blotting. Dot- and Slot-Immunoblotting. Immunohistochemistry. Other applications not tested. Optimal dilutions are dependent on conditions and should be determined by the user.
Specificity:	Reacts with the heavy and light chains of turkey IgG as demonstrated by ELISA. May also react with the light chains of other turkey immunoglobulins. Species: Turkey. Other species not tested.
Storage:	Store the antibody undiluted at 2-8°C. DO NOT FREEZE! This product is photosensitive and should be protected from light. Avoid repeated freezing and thawing. Shelf life: one year from despatch.

For research and in vitro use only. Not for diagnostic or therapeutic work.

Material Safety Datasheets are available at www.acris-antibodies.com or on request.

Antibody Hotline - Technical Questions - Antibody Location Service
Free Call: 0800-2274746 (Germany only) - www.acris-antibodies.com