

Polyclonal Antibody to GFP (Ads. to Hu, Ms, Rt Serum Proteins) - FITC

Alternate names: GFP-Tag, Green fluorescent protein

Catalog No.: AP09218FC-N

Quantity: 0.1 mg

Concentration: 1.0 mg/ml (by UV absorbance at 280 nm)

Background: Green fluorescence protein (GFP) is a 27 kDa protein derived from the jellyfish *Aequorea victoria*, which emits green light (emission peak at a wavelength of 509 nm) when excited by blue light (excitation peak at a wavelength of 395 nm). Green Fluorescent Protein (GFP) has become an invaluable tool in cell biology research, since its intrinsic fluorescence can be visualized in living cells. GFP fluorescence is stable under fixation conditions and suitable for a variety of applications. GFP has been widely used as a reporter for gene expression, enabling researchers to visualize and localize GFP-tagged proteins within living cells without the need for chemical staining. Other applications of GFP include assessment of protein-protein interactions through the yeast two hybrid system and measurement of distance between proteins through fluorescence energy transfer (FRET) protocols. GFP technology has considerably contributed to a greater understanding of cellular physiology.

YFP differs from GFP due to a mutation at T203Y; antibodies raised against full-length GFP should also detect YFP and other variants.

Uniprot ID: [P42212](#)

NCBI: [6100](#)

Host / Isotype: Rabbit / IgG

Immunogen: GST-Green Fluorescent Protein (GFP) fusion protein corresponding to the full length amino acid sequence (246aa) derived from the jellyfish *Aequorea victoria*.

Format: **State:** Lyophilized purified Ig fraction

Purification: Immunoaffinity Chromatography using Green Fluorescent Protein (*Aequorea victoria*) coupled to agarose beads followed by solid phase adsorption(s) to remove any unwanted reactivities.

Buffer System: 0.02M Potassium Phosphate, 0.15M Sodium Chloride, pH 7.2 with 10 mg/ml BSA (IgG and Protease free) as stabilizer and 0.01% (w/v) Sodium Azide as preservative.

Label: FITC – Fluorescein isothiocyanate (MW 390 daltons)

Absorption / Emission: 495 nm / 528 nm

Molar Ratio: 3.0 moles FITC per mole of Rabbit IgG

Reconstitution: Restore with 0.1 ml of deionized water (or equivalent).

For research and in vitro use only. Not for diagnostic or therapeutic work.

Material Safety Datasheets are available at www.acris-antibodies.com or on request.

Antibody Hotline - Technical Questions - Antibody Location Service
Free Call: 0800-2274746 (Germany only) - www.acris-antibodies.com

- Applications:** Suitable for Immunomicroscopy and Flow Cytometry or FACS analysis as well as other antibody based fluorescent assays requiring lot-to-lot consistency.
Recommended Dilutions:
Immunofluorescence: 1/500-1/2,500.
Flow Cytometry: 1/2,000-1/10,000.
Other applications not tested. Optimal dilutions are dependent on conditions and should be determined by the user.
- Specificity:** Assay by Immunoelectrophoresis resulted in a single precipitin arc against anti-Rabbit Serum, anti-Fluorescein and purified and partially purified Green Fluorescent Protein (*Aequorea victoria*) Serum.
No reaction was observed against Human, Mouse and Rat Serum Proteins.
- Storage:** Store vial at 2-8°C prior to restoration. Following restoration product can be stored undiluted at 2-8° for one month or (in aliquots) at -20°C or below.
Avoid repeated freezing and thawing. Centrifuge product if not completely clear after standing at room temperature.
Shelf life: One year from despatch.
- General References:** 1. The and Feltkamp, Immunology 18; 865, 1970.

For research and in vitro use only. Not for diagnostic or therapeutic work.

Material Safety Datasheets are available at www.acris-antibodies.com or on request.

Antibody Hotline - Technical Questions - Antibody Location Service
Free Call: 0800-2274746 (Germany only) - www.acris-antibodies.com