

Monoclonal Antibody to VZV / HHV-3 - FITC

Alternate names:	Human herpesvirus 3, Varicella-zoster virus
Catalog No.:	BM1346F
Quantity:	0.1 mg
Concentration:	1.3 mg/ml
Host / Isotype:	Mouse / IgM
Clone:	5B7
Immunogen:	Native, Ellen strain.
Format:	State: Liquid purified IgM Purification: Ammonium Sulphate fractionation Buffer System: PBS buffer, pH 7.2 Label: FITC
Applications:	Immunofluorescence: 0.1 - 2 µg/slide. Other applications not tested. Optimal dilutions are dependent on conditions and should be determined by the user.
Specificity:	Fluorescence is seen in both the nucleus and cytoplasm.
Storage:	The antibody can be shipped at 2-8°C. Store (in aliquots) at -20°C only. Avoid repeated freezing and thawing. Shelf life: one year from despatch.

For research and in vitro use only. Not for diagnostic or therapeutic work.

Material Safety Datasheets are available at www.acris-antibodies.com or on request.

Antibody Hotline - Technical Questions - Antibody Location Service
Free Call: 0800-2274746 (Germany only) - www.acris-antibodies.com