

Polyclonal Antibody to Human Plasminogen - FITC

Catalog No.:	BP749F
Quantity/Conc.:	1 ml / 5.1 mg/ml (OD 280 nm)
Immunogen:	Highly pure human Plasminogen
Host/Isotype:	Sheep
Format:	This antibody is supplied as liquid, FITC conjugated immunoglobulin fraction in Glycine buffered saline pH, 7.4, 0.09% NaN ₃ , 0.1% EACA, and 0.01% Benzamidine. Molar ratio FITC/protein is 2.1/1.0
Applications:	Immunofluorescence: 1/20 - 1/40 (as a guide) Other applications not tested. Optimal dilutions of this antibody are dependent on conditions and should be determined by the user.
Specificity:	Minimal reactivity with related serum proteins. May react with Plasminogen from other species.
Storage:	Aliquot and store the antibody at -20°C. Avoid repeated freezing and thawing. Protect from light. Shelf life: one year from despatch.

BP749F/CS0306

or research and *in vitro* use only. Not for diagnostic or therapeutic work.
Material Safety Datasheets are available at www.acris-antibodies.com or on request.

Antibody Hotline • Technical Questions • Antibody Location Service
Free Call: 0800-2274746 (Germany only)
service@acris-antibodies.de

www.acris-antibodies.com