

Monoclonal Antibody to Rat IgG2b (heavy chain) - Biotin

Catalog No.:	SM1497B
Quantity:	0.5 mg
Concentration:	1.0 mg/ml
Host / Isotype:	Mouse / IgG1
Clone:	MARG2b-8
Immunogen:	Rat IgG2b myeloma protein (IR863). Spleen cells from immunised BALB/c mice were fused with cells of the mouse SP2/0 myeloma cell line.
Format:	State: Liquid purified IgG Purification: Affinity chromatography Buffer System: PBS, pH7.4 containing 0.09% Sodium Azide and 50% Glycerol Label: Biotin
Applications:	ELISA: 500ng/ml. Other applications not tested. Optimal dilutions are dependent on conditions and should be determined by the user.
Specificity:	This antibody recognises IgG2b and shows no cross reaction with other immunoglobulin types. Species: Rat. Other species not tested.
Storage:	Store the antibody undiluted at 2-8°C for one month or (in aliquots) at -20°C for longer. Avoid repeated freezing and thawing. Shelf life: one year from despatch.

For research and in vitro use only. Not for diagnostic or therapeutic work.

Material Safety Datasheets are available at www.acris-antibodies.com or on request.

Antibody Hotline - Technical Questions - Antibody Location Service
Free Call: 0800-2274746 (Germany only) - www.acris-antibodies.com