

Monoclonal Antibody to Feline CD5 -FITC

Catalog No.:	SM1584F
Quantity/Conc.:	0.1 mg / 0.1 mg/ml
Clone:	FE1.1B11
Host/Isotype:	Mouse IgG1
Format:	This antibody is supplied as liquid, Protein G affinity purified immunoglobulin fraction, conjugated to Fluorescein Isothiocyanate Isomer 1 (FITC) in PBS with 0.09% sodium azide as preservative and 1% BSA as stabilizer.
Applications:	Flow cytometry: neat – 1/10; use 10 µl of the suggested working dilution to label 10 ⁶ cells. Other applications not tested. Optimal dilutions of this antibody are dependent on conditions and should be determined by the user.
Specificity:	This antibody recognises feline CD5, a 67kD cell surface glycoprotein expressed by peripheral T lymphocytes and thymocytes.
Storage:	Store the antibody at 4-8°C for one month or at –20°C for longer. This product is photosensitive and should be protected from light. This product should be stored undiluted. Avoid repeated freezing and thawing. Shelf life: one year from despatch.

SM1584F/ME0603

For research and *in vitro* use only. Not for diagnostic or therapeutic work.
Material Safety Datasheets are available at www.acris-antibodies.com or on request.

Antibody Hotline • Technical Questions • Antibody Location Service
0800-2274746
service@acris-antibodies.de

www.acris-antibodies.com