

Monoclonal Antibody to CD3 delta chain - FITC

Alternate names:	CD3D, T-cell surface glycoprotein CD3 delta chain, T3D
Catalog No.:	SM253F
Quantity:	0.1 mg
Concentration:	0.1 mg/ml
Background:	<p>T cell activation through the antigen receptor (TCR) involves the cytoplasmic tails of the CD3 subunits: CD3 gamma, CD3 delta, CD3 epsilon and CD3 zeta. These CD3 subunits are structurally related members of the immunoglobulins super family encoded by closely linked genes on human chromosome 11. The CD3 components have long cytoplasmic tails that associate with cytoplasmic signal transduction molecules. This association is mediated at least in part by a double tyrosine based motif present in a single copy in the CD3 subunits. CD3 may play a role in TCR induced growth arrest, cell survival and proliferation. The CD3 antigen is present on 68-82% of normal peripheral blood lymphocytes, 65-85% of thymocytes and Purkinje cells in the cerebellum. It is never expressed on B or NK cells. Decreased percentages of T lymphocytes may be observed in some autoimmune diseases.</p>
Uniprot ID:	P19377
NCBI:	NP_037301.1
GeneID:	25710
Host / Isotype:	Mouse / IgM
Clone:	1F4
Immunogen:	F344 rat T cells stimulated with PMA (TPA) and calcium ionophore.
Format:	<p>State: Liquid purified IgM fraction. Buffer System: PBS, pH 7.4 containing 0.09% Sodium Azide as preservative and 1% BSA as stabilizer. Label: FITC – Fluorescein Isothiocyanate Isomer 1</p>
Applications:	<p>Flow cytometry: Use 10 µl of Neat antibody to label 10e6 cells in 100 µl. Functionally the addition of the antibody to a culture of rat T cells induces the proliferation of T-cells in the presence of PMA (2). We recommend antibody SM253A for this purpose. Other applications not tested. Optimal dilutions are dependent on conditions and should be determined by the user.</p>
Specificity:	<p>SM253F recognises rat CD3, a 25 kD antigen which is found on rat T cells. SM253F does not react with rat B cells. Species: Rat. Other species not tested.</p>

For research and in vitro use only. Not for diagnostic or therapeutic work.

Material Safety Datasheets are available at www.acris-antibodies.com or on request.

Antibody Hotline - Technical Questions - Antibody Location Service
Free Call: 0800-2274746 (Germany only) - www.acris-antibodies.com

Storage:

Store the antibody undiluted at 2-8°C or (in aliquots) at -20°C for longer.
This product is photosensitive and should be protected from light.
Avoid repeated freezing and thawing.
Shelf life: one year from despatch.

General References:

1. Nicolls, M.G. et al. (1992). Induction of long term specific tolerance to allografts in vivo by therapy with a monoclonal antibody to rat CD3. *Transplantation* 55: 459-468.
2. Tanaka, T. et al. (1989). Characterization of a CD3-like rat T cells surface antigen recognised by a monoclonal antibody. *J. Immunol.* 142: 2791-2795.

For research and in vitro use only. Not for diagnostic or therapeutic work.

Material Safety Datasheets are available at www.acris-antibodies.com or on request.

Antibody Hotline - Technical Questions - Antibody Location Service
Free Call: 0800-2274746 (Germany only) - www.acris-antibodies.com