

Monoclonal Antibody to Porcine CD8 -FITC

Catalog No.:	SM486F
Quantity/Conc.:	0.1 mg / 0.1 mg/ml
Clone:	MIL-12
Host/Isotype:	Mouse IgG2a
Immunogen:	Porcine mesenteric lymphocytes.
Format:	This antibody is supplied as liquid, Protein G affinity purified immunoglobulin fraction, conjugated to Fluorescein Isothiocyanate Isomer 1 (FITC) in PBS buffer pH 7.4 with 0.1% sodium azide as preservative and 1% BSA as stabilizer. FITC:Protein (molar ratio): 6.3 : 1.0.
Applications:	Flow cytometry: use 10 µl of neat antibody to label 10 ⁶ cells. Other applications not tested.
Specificity:	This antibody recognises a subset of porcine T-cells. Cross reactivities not tested.
Storage:	Store the antibody at 4-8°C for one month or at -20°C for longer. This product is photosensitive and should be protected from light. Should this product contain a precipitate we recommend microcentrifugation before use. Avoid repeated freezing and thawing. Shelf life: one year from despatch.

SM486F/ME1102

For research and *in vitro* use only. Not for diagnostic or therapeutic work.
Material Safety Datasheets are available at www.acris-antibodies.com or on request.

Antibody Hotline • Technical Questions • Antibody Location Service
0800-2274746
service@acris-antibodies.de

www.acris-antibodies.com