

Monoclonal Antibody to Porcine CD8 - PE

Catalog No.:	SM486R
Quantity/Conc.:	100 Tests (lyophil.)
Clone:	MIL-12
Host/Isotype:	Mouse IgG2a
Immunogen:	Porcine mesenteric lymphocytes.
Format:	This antibody is supplied lyophilized from a Protein G affinity purified immunoglobulin fraction, conjugated to R. Phycoerythrin (RPE) in PBS buffer pH 7.4 with 0.09% sodium azide as preservative and 1% BSA as stabilizer. Reconstitute with 1.0 ml distilled water
Applications:	Flow cytometry: use 10 µl of neat antibody to label 10 ⁶ cells in 100 µl. Other applications not tested. Optimal dilutions are dependent on conditions and should be determined by the user.
Specificity:	This antibody recognises a subset of porcine T-cells. Cross reactivities not tested.
Storage:	Store lyophilized and undiluted reconstituted antibody at +4°C. DO NOT FREEZE. This product is photosensitive and should be protected from light. Should this product contain a precipitate we recommend microcentrifugation before use. Shelf life: one year from despatch.

SM486R/PB0904

For research and *in vitro* use only. Not for diagnostic or therapeutic work.
Material Safety Datasheets are available at www.acris-antibodies.com or on request.

Antibody Hotline • Technical Questions • Antibody Location Service
0800-2274746
service@acris-antibodies.de

www.acris-antibodies.com