



Acris Antibodies GmbH
Im Himmelreich 11
32120 Hiddenhausen
Phone: +49-5221-34606-0
Fax: +49-5221-34606-11
info@acris-online.de

Monoclonal Antibody to Pig CD45RA -FITC

Catalog No.:	SM563F
Quantity/Conc.:	0.1 mg / 0.1 mg/ml
Clone:	MIL13
Host/Isotype:	Mouse IgG1
Immunogen:	Cells isolated from porcine mesenteric lymph node
Format:	This antibody is supplied as liquid, purified immunoglobulin fraction, conjugated to Fluorescein Isothiocyanate Isomer 1 (FITC) with 0.09% sodium azide as preservative and 1% BSA as stabilizer.
Applications:	Flow cytometry, use 10 µl of neat antibody to label 10 ⁶ cells. Optimal dilutions of this antibody are dependent on application and should be determined by the user.
Specificity:	MIL13 recognises the porcine CD45RA cell surface antigen. Species Cross Reactivity: None known
Storage:	Store the antibody at 4-8°C for one month or at -20°C for longer. This product is photosensitive and should be protected from light. Avoid repeated freezing and thawing. Shelf life: one year from despatch.

SM563F/ME0803

For research and *in vitro* use only. Not for diagnostic or therapeutic work.
Material Safety Datasheets are available at www.acris-antibodies.com or on request.

Antibody Hotline - Technical Questions - Antibody Location Service
0800-2274746
service@acris-antibodies.de

www.acris-antibodies.com