



Acris Antibodies GmbH  
Im Himmelreich 11  
32120 Hiddenhausen  
Phone: +49-5221-34606-0  
Fax: +49-5221-34606-11  
info@acris-online.de

---

## Monoclonal Antibody to Pig CD45RC

<b>Catalog No.:</b>	SM564
<b>Quantity/Conc.:</b>	2.0 ml
<b>Clone:</b>	MIL15
<b>Host/Isotype:</b>	Mouse IgM
<b>Immunogen:</b>	Cells isolated from porcine mesenteric lymph node
<b>Format:</b>	This antibody is supplied as liquid tissue culture supernatant with 0.1% sodium azide as preservative.
<b>Applications:</b>	Flow cytometry, use 10 µl of neat antibody to label 10 <sup>6</sup> cells. Immunohistochemistry: frozen sections and paraffin sections. Optimal dilutions of this antibody are dependent on application and should be determined by the user.
<b>Specificity:</b>	MIL15 recognises the porcine CD45RC cell surface antigen. Species Cross Reactivity: None known
<b>Storage:</b>	Store the antibody at 4-8°C for one month or at -20°C for longer. Avoid repeated freezing and thawing. Shelf life: one year from despatch.

SM564/DJS1101

**For research and *in vitro* use only. Not for diagnostic or therapeutic work.**  
**Material Safety Datasheets are available at [www.acris-antibodies.com](http://www.acris-antibodies.com) or on request.**

---

Antibody Hotline • Technical Questions • Antibody Location Service  
0800-2274746  
service@acris-antibodies.de

[www.acris-antibodies.com](http://www.acris-antibodies.com)