

Goat Polyclonal Antibody to Rat IgG (Mouse Absorbed)

Catalog No.:	SP1020P
Quantity:	0.5 mg
Concentration:	0.5 mg/ml
Host:	Goat
Immunogen:	Rat IgG
Format:	State: Liquid purified IgG Purification: Affinity chromatography Buffer System: PBS containing 0.09% Sodium Azide
Applications:	Flow Cytometry: 5µg/ml. Immunoprecipitation: 1µg/ml - 5µg/ml. Other applications not tested. Optimal dilutions are dependent on conditions and should be determined by the user.
Specificity:	This antibody recognises IgG. Cross-reactivity with mouse IgG has been minimised by adsorption. Species: Rat. Other species not tested.
Storage:	Store the antibody undiluted at 2-8°C for one month or (in aliquots) at -20°C for longer. Avoid repeated freezing and thawing. Shelf life: one year from despatch.

For research and in vitro use only. Not for diagnostic or therapeutic work.

Material Safety Datasheets are available at www.acris-antibodies.com or on request.

Antibody Hotline - Technical Questions - Antibody Location Service
Free Call: 0800-2274746 (Germany only) - www.acris-antibodies.com