

Fab Fragment of Rabbit anti-Rat IgG

Catalog No.:	SP1021P
Quantity:	1 mg
Concentration:	1.0 mg/ml
Host:	Rabbit
Immunogen:	Purified polyclonal rat IgG
Format:	State: Liquid Fab fragment of IgG Purification: Affinity chromatography Buffer System: PBS, pH 7.4 containing 0.09% Sodium Azide
Applications:	ELISA: 1/250 - 1/1000. Immunodiffusion. Other applications not tested. Optimal dilutions are dependent on conditions and should be determined by the user.
Specificity:	Fab fragment recognises rat IgG, detecting all subclasses in immunodiffusion assays. Species: Rat. Other species not tested.
Storage:	Store the antibody undiluted at 2-8°C for one month or (in aliquots) at -20°C for longer. Avoid repeated freezing and thawing. Shelf life: one year from despatch.

For research and in vitro use only. Not for diagnostic or therapeutic work.

Material Safety Datasheets are available at www.acris-antibodies.com or on request.

Antibody Hotline - Technical Questions - Antibody Location Service
Free Call: 0800-2274746 (Germany only) - www.acris-antibodies.com