

Polyclonal Antibody to Bovine IgG1 - FITC

Alternate names:	Immunoglobulin G1
Catalog No.:	SP1049F
Quantity:	1 mg
Concentration:	1.0 mg/ml
Host:	Sheep
Immunogen:	Purified bovine IgG1
Format:	State: Liquid IgG fraction Purification: Affinity chromatography Buffer System: PBS containing 0.09% Sodium Azide Label: FITC – Fluorescein Isothiocyanate Isomer 1
Applications:	Flow Cytometry: 1/20 - 1/100. Immunohistochemistry on frozen sections: 1/20 - 1/100. Other applications not tested. Optimal dilutions are dependent on conditions and should be determined by the user.
Specificity:	This antibody recognises IgG1 and shows no cross - reactivity with other immunoglobulin classes in immunoelectrophoresis. Species: Bovine. Other species not tested.
Storage:	Store the antibody undiluted at 2-8°C for one month or (in aliquots) at -20°C for longer. Avoid repeated freezing and thawing. This product is photosensitive and should be protected from light. Shelf life: one year from despatch.

For research and in vitro use only. Not for diagnostic or therapeutic work.

Material Safety Datasheets are available at www.acris-antibodies.com or on request.

Antibody Hotline - Technical Questions - Antibody Location Service
Free Call: 0800-2274746 (Germany only) - www.acris-antibodies.com